

Saneso™

	Diameter (mm)	Instrument channel (mm)	Working length (mm)	Total length (mm)	Field of view (°)	Observation range (mm)	Bending angle vertical/horizontal	Display mode	Select Band Imaging (TM)	Dual Band Imaging (TM)
Colonoscope										
SAN 360C-A: Saneso 360 Colonoscope Adult	13.8	3.7	1680	3010	360°	2 - 100	180°/160°	Panoramic integrated	y	y
SAN 1C-A: Saneso Colonoscope Slim	11.8	3.7	1680	3010	160°	2 - 100	180°/160°	Traditional forward only	y	n

Gastroscope										
SAN 360G-A: Saneso 360 Gastroscope	13.8	3.7	1030	2360	360°	2 - 100	180°/160°	Panoramic integrated	y	y
SAN 1G-A: Saneso Gastroscope	11.8	3.7	1030	2360	160°	2 - 100	180°/160°	Traditional forward only	y	n

Feature	Saneso™	Traditional
Revolutionary		
Maximum FOV - colonoscope	360°	140° -170°
Maximum FOV - gastroscope	360°	140°
Image	Single panoramic	Single limited
Spectrum / bands	White and Narrow Spectra	White and NBI
Simultaneous multi-band imaging	Dual Band Imaging™	None
Lesion measurement assist	Calibration Markers	None
Lens wash	5	1
Ease of Use	Simple	Simple
Ergonomic improvement	Yes (comfi-grip™)	-

Familiar operation		
Instrument channel (colonoscope)	3.7 mm	3.2 / 3.7 mm
Instrument channel (gastroscope)	3.7 mm	2.8 / 3.2 / 3.7 mm
Flexibility colonoscope (V/H)	180° / 160°	180° / 160°
Flexibility gastroscope (V/H)		210° up & 90° down /100°
Stiffness	Selective Stiffness™	Variable Stiffness
Diameter (colonoscope)	13.8 / 11.8 mm	13.2 / 11.7 mm
Diameter (gastroscope)	13.8 / 11.8 mm	12.6 / 9.9 mm
Dedicated tissue wash channel	Yes	Yes
Distal head length	3 cm	3 cm
Image	HD	HD
Display	ultra HD	HD
Depth of field	2-100 mm	2-100 mm

Affordable		
Cost	CALL	***
Warranty	24 months	12 months
Peace of mind	360 extended coverage (4 years)	Variable
60 day money back	Yes	No
Lifetime cost of ownership	**	***